



Alasta™

ALOE

Soft Nitrile Meets Aloe Vera to Create Alasta Aloe

- The nitrile glove with aloe designed for latex users
- Soft Fit™ nitrile gloves
- Aloe vera coated inner surface
- Performs like latex
- Economically priced
- Extreme comfort with precision fit
- Textured fingers

Available in X-Small, Small, Medium, Large, X-Large

DASH®
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Alasta™ Aloe Glove Features

Material: acrylonitrile butadiene polymer
Type: non-sterile
Color: green
Glove Interior: powder free, aloe vera coated
Glove Exterior: textured fingers
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
 Before Aging (MPa): 14
 After Aging (MPa): 14
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	74 (±4)	250 (±10)	0.08 mm (±.015)/3.2 mil	0.10 mm (±.02)/3.9 mil
Small	81 (±4)	250 (±10)	0.08 mm (±.015)/3.2 mil	0.10 mm (±.02)/3.9 mil
Medium	97 (±4)	250 (±10)	0.08 mm (±.015)/3.2 mil	0.10 mm (±.02)/3.9 mil
Large	107 (±4)	250 (±10)	0.08 mm (±.015)/3.2 mil	0.10 mm (±.02)/3.9 mil
X-Large	115 (±4)	250 (±10)	0.08 mm (±.015)/3.2 mil	0.10 mm (±.02)/3.9 mil