

AloePRO™ Synthetic



Soothe Your Hands and Ease Your Mind.

DASH AloePRO Synthetic exam gloves are your latex-free solution to better feeling gloves.

- Active Aloe™ inner surface soothes dry irritated hands by improving your skin's ability to hydrate itself
- Made from innovative butadiene technology with stress relaxation properties that are superior to nitrile
- Priced lower than most nitrile gloves
- Far superior fit than vinyl gloves
- Six sizes for a precise fit

Available in X-Small, Small, Small-Medium, Medium, Large, X-Large

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&T Office

AloePRO™ Synthetic Glove Features

Material: butadiene
Type: non-sterile
Color: green
Glove Interior: powder free
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9001, FDA QSR
Tensile Strength
 Before Aging (MPa): 14
 After Aging (MPa): 14
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	74 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Small	81 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Small-Medium	89 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Medium	97 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Large	107 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil
X-Large	115 (±4)	250 (±10)	0.10 mm (±0.015)/3.9 mil	0.125 mm (±0.02)/4.9 mil