

Black Maxx™ Nitrile



Protection to the Maxx.

- Deep-Black color with a unique gloss finish
- Fully-textured surface provides a secure grip even when hands are wet
- New technologically advanced nitrile formulation creates a softer and more elastic glove while maintaining strength and durability
- Softer fit decreases hand fatigue, increases dexterity and provides a greater range of motion

Available in Small, Medium, Large, X-Large, XX-Large

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Black Maxx™ Nitrile Glove Features

Material: acrylonitrile butadiene polymer
Type: non-sterile
Color: black
Glove Interior: powder free
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9001, ISO 13485
Tensile Strength
 Before Aging (MPa): 14 min
 After Aging (MPa): 14 min
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
Small	85 (±5)	240	0.1 mm (±0.03)/3.9 mil	0.12 mm (±0.03)/4.7 mil
Medium	95 (±5)	240	0.1 mm (±0.03)/3.9 mil	0.12 mm (±0.03)/4.7 mil
Large	106 (±5)	240	0.1 mm (±0.03)/3.9 mil	0.12 mm (±0.03)/4.7 mil
X-Large	115 (±5)	240	0.1 mm (±0.03)/3.9 mil	0.12 mm (±0.03)/4.7 mil
XX-Large	126 (±5)	240	0.1 mm (±0.03)/3.9 mil	0.12 mm (±0.03)/4.7 mil