


Black Maxx™



Protection to the Maxx.

- Jet black gloves with a unique gloss finish
- New micro-textured surface provides a secure grip even when hands are wet
- Premium latex offers tactile sensitivity and elasticity for a soft, comfortable fit over long periods of wear
- Impressive tear resistance and strength created by slightly thicker gloves helps maintain their durability

Available in X-Small, Small, Medium, Large, X-Large

DASH®
MEDICAL GLOVES

Where Safety Meets Comfort™

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Black Maxx™ Glove Features

Material: natural rubber latex
Type: non-sterile
Color: black
Glove Interior: powder free
Glove Exterior: textured fingers
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: less than 50 µg/g per ASTM 5712
Quality Certification: ISO 9001, ISO 13485
Tensile Strength
 Before Aging (MPa): 18 min
 After Aging (MPa): 14 min
Elasticity/Elongation
 Before Aging: 650%
 After Aging: 500%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	75 (±5)	240	0.14 mm (±0.03)/5.5 mil	0.16 mm (±0.03)/6.3 mil
Small	85 (±5)	240	0.14 mm (±0.03)/5.5 mil	0.16 mm (±0.03)/6.3 mil
Medium	95 (±5)	240	0.14 mm (±0.03)/5.5 mil	0.16 mm (±0.03)/6.3 mil
Large	106 (±5)	240	0.14 mm (±0.03)/5.5 mil	0.16 mm (±0.03)/6.3 mil
X-Large	115 (±5)	240	0.14 mm (±0.03)/5.5 mil	0.16 mm (±0.03)/6.3 mil