

FUTURA™

In Harmony With Your Skin's Natural Physiology™

Introducing DASH Futura — formulated to match your hand's natural pH level of 5.5

Products that come in contact with your skin often disrupt its natural mildly acidic environment. DASH Futura gloves, however, work with your skin to maintain a healthy pH level of 5.5, thereby preserving your body's first line of defense against irritants and infection. That's why DASH Futura gloves help relieve the stress on your skin caused by overwashing and frequent gloving.

Available in X-Small, Small, Medium, Large, X-Large



DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

FUTURA™ Glove Features

Material: natural rubber latex
Type: non-sterile
Color: teal blue
Glove Interior: powder free, pH balance formulated
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: less than 50 µg/g per ASTM D 5712-95
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
Before Aging (MPa): 24
After Aging (MPa): 20
Elasticity/Elongation
Before Aging: 750%
After Aging: 700%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	75 (±3)	250 (±10)	0.14 mm (±0.025)/5.6 mil	0.155 mm (±0.03)/6.1 mil
Small	85 (±3)	250 (±10)	0.14 mm (±0.025)/5.6 mil	0.155 mm (±0.03)/6.1 mil
Medium	96 (±3)	250 (±10)	0.14 mm (±0.025)/5.6 mil	0.155 mm (±0.03)/6.1 mil
Large	107 (±3)	250 (±10)	0.14 mm (±0.025)/5.6 mil	0.155 mm (±0.03)/6.1 mil
X-Large	117 (±3)	250 (±10)	0.14 mm (±0.025)/5.6 mil	0.155 mm (±0.03)/6.1 mil