

halo™



A Better Glove...By Design.

DASH halo gloves are the best performing nitrile gloves you can buy. halo uses the newest, most advanced nitrile formulation, offering superior tactile sensitivity for the precision your job requires.

- Super soft, rubber-like feel
- Durable, puncture resistant and protein free
- Consistent high quality
- Micro-textured fingers for enhanced grip

Available in X-Small, Small, Medium, Large, X-Large

DASH®
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&T Office

halo™ Glove Features

Material: acrylonitrile butadiene polymer
Type: non-sterile
Color: blue
Glove Interior: powder free
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
 Before Aging (MPa): 18
 After Aging (MPa): 16
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	76 (±4)	240	0.09 mm (±0.02)/3.5 mil	0.12 mm (±0.02)/4.7 mil
Small	86 (±4)	240	0.09 mm (±0.02)/3.5 mil	0.12 mm (±0.02)/4.7 mil
Medium	98 (±4)	240	0.09 mm (±0.02)/3.5 mil	0.12 mm (±0.02)/4.7 mil
Large	107 (±4)	240	0.09 mm (±0.02)/3.5 mil	0.12 mm (±0.02)/4.7 mil
X-Large	115 (±4)	240	0.09 mm (±0.02)/3.5 mil	0.12 mm (±0.02)/4.7 mil