

Hi-Risk Protector™ Latex



Premium Protection for High-Risk Conditions.

Extra cuff length, extra thickness, durable latex—all make DASH HRP Latex exam gloves premium protection in high-risk conditions.

- Thicker glove
- Extended cuff
- Durable latex
- Powder free
- Fully textured grip

Available in Small, Medium, Large, X-Large, XX-Large

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort[™]

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Hi-Risk Protector™ Latex Glove Features

Material: natural rubber latex
Type: non-sterile
Color: natural
Glove Interior: powder free, polymer coated
Glove Exterior: textured fingers
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 50 gloves per dispenser box, by weight

Protein Content: less than 50 µg/g per ASTM D 3578-05
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
 Before Aging (MPa): 18 min
 After Aging (MPa): 14 min
Elasticity/Elongation
 Before Aging: 650%
 After Aging: 500%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
Small	85 (±5)	295	0.23 mm/9.1 mil	0.36 mm/14 mil
Medium	95 (±5)	295	0.23 mm/9.1 mil	0.36 mm/14 mil
Large	106 (±5)	295	0.23 mm/9.1 mil	0.36 mm/14 mil
X-Large	116 (±5)	295	0.23 mm/9.1 mil	0.36 mm/14 mil
XX-Large	118 (±5)	295	0.23 mm/9.1 mil	0.36 mm/14 mil