

Hi-Risk Protector™



Premium Protection for High-Risk Conditions.

Extra cuff length, extra thickness, puncture-resistant nitrile—all make DASH HRP exam gloves premium protection in high-risk conditions.

- Thicker glove
- Extended cuff
- Puncture resistant nitrile
- Latex free and powder free
- Textured fingers

Available in Small, Medium, Large, X-Large, XX-Large

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Hi-Risk Protector™ Glove Features

Material: acrylonitrile butadiene polymer
Type: non-sterile
Color: steel blue
Glove Interior: powder free
Glove Exterior: textured fingers
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 50 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
 Before Aging (MPa): 14
 After Aging (MPa): 14
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
Small	85 (±5)	295 (±5)	0.14 mm (±0.015)/5.5 mil	0.18 mm (±0.015)/7.1 mil
Medium	95 (±5)	295 (±5)	0.14 mm (±0.015)/5.5 mil	0.18 mm (±0.015)/7.1 mil
Large	106 (±5)	295 (±5)	0.14 mm (±0.015)/5.5 mil	0.18 mm (±0.015)/7.1 mil
X-Large	113 (±5)	295 (±5)	0.14 mm (±0.015)/5.5 mil	0.18 mm (±0.015)/7.1 mil
XX-Large	120 (±5)	295 (±5)	0.14 mm (±0.015)/5.5 mil	0.18 mm (±0.015)/7.1 mil