

LP™

Fit for a King. Priced for Much Less.

Take comfort in knowing there are exam gloves perfectly sized for your hands. Available in six sizes, DASH LP is the professional's solution for exceptional comfort and fit. DASH LP exam gloves are made from high quality, low modulus latex which provides a barrier that is so form fitting, you hardly notice you're wearing gloves at all.

Available in X-Small (6), Small (6½), Small-Medium (7), Medium (7½), Medium-Large (8), Large (8½)



DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™

CAUTION: This product contains natural rubber latex which may cause allergic reactions.



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

LP™ Glove Features

Material: natural rubber latex
Type: non-sterile
Color: natural
Glove Interior: powdered with USP grade Absorbo HP™
Glove Exterior: smooth surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: less than 200 µg/g per ASTM D 5712-99
Quality Certification: ISO 9002, FDA QSR
Tensile Strength
Before Aging (MPa): 21 min
After Aging (MPa): 16 min
Elasticity/Elongation
Before Aging: 700%
After Aging: 500%

tGlove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	76 (±4)	240 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil
Small	84 (±4)	240 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil
Small-Medium	88 (±4)	245 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil
Medium	94 (±4)	245 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil
Medium-Large	102 (±4)	245 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil
Large	110 (±4)	245 (±5)	0.13 mm (±0.04)/5.1 mil	0.14 mm (±0.04)/5.5 mil