

## Nitrile PF with Aloe™



## Soothingly Tough From the Inside Out.

Puncture resistant and fully textured for an enhanced grip, DASH Nitrile PF with Aloe gloves are the perfect solution for dry, irritated hands that perform high-risk jobs. As durable nitrile protects your hands from outside contact, the Active Aloe™ inner surface soothes hands by holding in your skin's moisture. Protein free and powder free, DASH Nitrile PF with Aloe gloves are also tested and approved for use with chemotherapy drugs.

Available in X-Small, Small, Medium, Large, X-Large

**DASH**<sup>®</sup>  
MEDICAL GLOVES

*Where Safety Meets Comfort™*



Phone: 1-800-523-2055 Fax: 1-800-523-7795 [www.dashmedical.com](http://www.dashmedical.com)

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

## Nitrile PF with Aloe™ Glove Features

**Material:** acrylonitrile butadiene polymer  
**Type:** non-sterile  
**Color:** green  
**Glove Interior:** powder free, aloe vera coated  
**Glove Exterior:** textured surface  
**Shape:** ambidextrous  
**Cuff:** beaded  
**Usage:** single use only  
**Quantities:** 100 gloves per dispenser box, by weight

**Protein Content:** none  
**Quality Certification:** ISO 9002, FDA QSR  
**Tensile Strength**  
**Before Aging (MPa):** 14 min  
**After Aging (MPa):** 14 min  
**Elasticity/Elongation**  
**Before Aging:** 500%  
**After Aging:** 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	74 (±4)	250 (±10)	0.125 mm (±0.02)/4.9 mil	0.15 mm (±0.02)/5.9 mil
Small	81 (±4)	250 (±10)	0.125 mm (±0.02)/4.9 mil	0.15 mm (±0.02)/5.9 mil
Medium	96 (±4)	250 (±10)	0.125 mm (±0.02)/4.9 mil	0.15 mm (±0.02)/5.9 mil
Large	106 (±4)	250 (±10)	0.125 mm (±0.02)/4.9 mil	0.15 mm (±0.02)/5.9 mil
X-Large	115 (±4)	250 (±10)	0.125 mm (±0.02)/4.9 mil	0.15 mm (±0.02)/5.9 mil