

blue pearlescent Nitrile PF™



Protection with *Style!*

Go ahead—brighten your workday with the shimmer of blue pearlescent Nitrile PF gloves. Much more than dazzling looks, these gloves are made from the latest generation of soft nitrile material to provide the premium protection nitrile is known for, while avoiding the restrictiveness associated with traditional nitrile gloves.

Comfort, protection and style. You can't go wrong when choosing DASH blue pearlescent Nitrile PF gloves for your hand protection.

Available in X-Small, Small, Medium, Large, X-Large, XX-Large



DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort[™]



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2009 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

blue pearlescent Nitrile PF™ Glove Features

Material: acrylonitrile butadiene polymer
Type: non-sterile
Color: pearlescent blue
Glove Interior: powder free
Glove Exterior: textured surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9001, ISO 13485, cGMP, CE
Tensile Strength
 Before Aging (MPa): 14 min
 After Aging (MPa): 14 min
Elasticity/Elongation
 Before Aging: 500%
 After Aging: 400%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	76 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Small	86 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Medium	97 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil
Large	107 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil
X-Large	115 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil
XX-Large	124 (±4)	250 (±10)	0.10 mm (±0.02)/3.9 mil	0.125 mm (±0.02)/4.9 mil