

Nuvo™

So Soft and Flexible You Won't Believe It's Vinyl.

These are, hands down, the most amazing vinyl exam gloves you'll find. Slide on a pair of DASH Nuvo white vinyl gloves and see for yourself how soft and form fitting they feel. DASH Nuvo gloves are specially designed to retain the superior elongation properties found in natural latex yet normally absent from most synthetics such as standard vinyl.

Available in X-Small, Small, Medium, Large, X-Large



DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Nuvo™ Glove Features

Material: vinyl
Type: non-sterile
Color: white
Glove Interior: powder free
Glove Exterior: smooth surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9001, ISO 13485, FDA QSR
Tensile Strength
 Before Aging (MPa): 11 min
 After Aging (MPa): 11 min
Elasticity/Elongation
 Before Aging: 350%
 After Aging: 350%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	75 (±5)	240 (±5)	0.11 mm (±0.03)/4.3 mil	0.09 mm (±0.03)/3.5 mil
Small	85 (±5)	240 (±5)	0.11 mm (±0.03)/4.3 mil	0.09 mm (±0.03)/3.5 mil
Medium	95 (±5)	240 (±5)	0.11 mm (±0.03)/4.3 mil	0.09 mm (±0.03)/3.5 mil
Large	105 (±5)	240 (±5)	0.11 mm (±0.03)/4.3 mil	0.09 mm (±0.03)/3.5 mil
X-Large	115 (±5)	240 (±5)	0.11 mm (±0.03)/4.3 mil	0.09 mm (±0.03)/3.5 mil