

# Vinyl LP™



## Versatility at a Comfortable Price.

If you need a versatile alternative to latex gloves at a comfortable price, DASH Vinyl LP exam gloves are your answer. Latex-free DASH Vinyl gloves are the most economical way to protect your hands while eliminating concerns of latex protein allergies. DASH Vinyl gloves also work great when mixing impressions or other materials that may stick to latex gloves.

Available in X-Small, Small, Medium, Large, X-Large

**DASH**<sup>®</sup>  
MEDICAL GLOVES

*Where Safety Meets Comfort™*



Phone: 1-800-523-2055 Fax: 1-800-523-7795 [www.dashmedical.com](http://www.dashmedical.com)

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

## Vinyl LP™ Glove Features

**Material:** vinyl  
**Type:** non-sterile  
**Color:** clear  
**Glove Interior:** powdered with USP grade Absorbo HP™  
**Glove Exterior:** smooth surface  
**Shape:** ambidextrous  
**Cuff:** beaded  
**Usage:** single use only  
**Quantities:** 100 gloves per dispenser box, by weight

**Protein Content:** none  
**Quality Certification:** ISO 13485, FDA QSR  
**Tensile Strength**  
**Before Aging (MPa):** 11 min  
**After Aging (MPa):** 11 min  
**Elasticity/Elongation**  
**Before Aging:** 300%  
**After Aging:** 300%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	80 (±5)	230 (±10)	0.08 mm (±0.02)/3.2 mil	0.08 mm (±0.02)/3.2 mil
Small	85 (±5)	230 (±10)	0.08 mm (±0.02)/3.2 mil	0.08 mm (±0.02)/3.2 mil
Medium	95 (±5)	230 (±10)	0.08 mm (±0.02)/3.2 mil	0.08 mm (±0.02)/3.2 mil
Large	105 (±5)	230 (±10)	0.08 mm (±0.02)/3.2 mil	0.08 mm (±0.02)/3.2 mil
X-Large	115 (±5)	230 (±10)	0.08 mm (±0.02)/3.2 mil	0.08 mm (±0.02)/3.2 mil