

Vinyl PF™



Versatility at a Comfortable Price.

If you need a versatile alternative to latex gloves at a comfortable price, DASH Vinyl PF exam gloves are your answer. Latex-free DASH Vinyl gloves are the most economical way to protect your hands while eliminating concerns of latex protein allergies. DASH Vinyl gloves also work great when mixing impressions or other materials that may stick to latex gloves.

Available in X-Small, Small, Medium, Large, X-Large

DASH[®]
MEDICAL GLOVES

Where Safety Meets Comfort™



Phone: 1-800-523-2055 Fax: 1-800-523-7795 www.dashmedical.com

©2008 DASH Medical Gloves Inc. DASH is a registered trademark of DASH Medical Gloves Inc., registered in the US P&TM Office

Vinyl PF™ Glove Features

Material: vinyl
Type: non-sterile
Color: clear
Glove Interior: powder free
Glove Exterior: smooth surface
Shape: ambidextrous
Cuff: beaded
Usage: single use only
Quantities: 100 gloves per dispenser box, by weight

Protein Content: none
Quality Certification: ISO 9001, ISO 13485, FDA QSR
Tensile Strength
 Before Aging (MPa): 11 min
 After Aging (MPa): 11 min
Elasticity/Elongation
 Before Aging: 400%
 After Aging: 350%

Glove Size	Palm Width (mm)	Length (mm)	Thickness at Palm	Thickness at Finger
X-Small	75 (±5)	240 (±5)	0.11 mm (±0.02)/4.3 mil	0.09 mm (±0.02)/3.5 mil
Small	85 (±5)	245 (±5)	0.11 mm (±0.02)/4.3 mil	0.09 mm (±0.02)/3.5 mil
Medium	95 (±5)	245 (±5)	0.11 mm (±0.02)/4.3 mil	0.09 mm (±0.02)/3.5 mil
Large	105 (±5)	245 (±5)	0.11 mm (±0.02)/4.3 mil	0.09 mm (±0.02)/3.5 mil
X-Large	115 (±5)	250 (±5)	0.11 mm (±0.02)/4.3 mil	0.09 mm (±0.02)/3.5 mil