

BlackMaxx®

Latex exam gloves



Comfort of Latex in a Bold Color



- Deep black color masks stains
- Excellent range of motion while
- Strong and durable
- Powder free
- Fully textured



"LOVEEEE"

Our staff LOVE these gloves! Easy, No-slip grip for our tattooing and piercing!"

- DASH Website Review



Bold black color provides a framing effect for better color definition and camouflaging stains on the glove eliminating distractions. Fully textured surface and advanced double coating techniques provides a secure grip even when wet.



PRODUCT DETAILS	
Material	Latex
Interior	Powder Free
Grade	Medical Examination
Watertight Test	1.5 AQL
Color	Black
Texture	Full

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.14 mm 5.5 mil
Finger Thickness	0.17 mm 6.6 mil
Length	240 mm 9.5 in
Tensile Strength (MPA)	18
Elongation %	650

ORDER INFORMATION	
Product Number	BMX100
Available Sizes	XS, S, M, L, XL
Glove Count	100/box 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested [†]	No
Tested Against Fentanyl	No
PFAS Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D3578, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Dental, Healthcare, Tattoo, Industrial Safety, Automotive

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings



*Typical values have an allowable specification variance
[†]Contact DASH Medical for a list of tested Chemotherapy drugs