

# GloveUp<sup>®</sup> 300

Nitrile exam gloves



## Flexible Thin Formulation, Value Packed



- Thin, flexible nitrile formulation
- Cobalt blue color
- Value packed 300 gloves per box
- Powder free
- Textured fingers
- Tested against fentanyl



### "Very Soft, Fits My Hand.

This glove is very soft, fits my hand like a real GLOVE. We also love the texture on the fingers. It now comes in a 300 count box, which saves us a lot since we use gloves like tissue paper!"

– ALEX



**Super light formulation** conforms to your hand effortlessly and delivers superior tactile sensitivity.



PRODUCT DETAILS	
Material	Nitrile
Interior	Powder Free
Grade	Medical Examination
Watertight Test	1.5 AQL
Color	Cobalt
Texture	Finger

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.05 mm   1.9 mil
Finger Thickness	0.08 mm   3.1 mil
Length	230 mm   9.0 in
Tensile Strength (MPA)	18
Elongation %	500

ORDER INFORMATION	
Product Number	GUP300
Available Sizes	XS, S, M, L, XL
Glove Count	300/box (XL 250/box) 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested <sup>†</sup>	Yes
Tested Against Fentanyl	Yes
PFAS Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D6319, ASTM D6978 (Chemo and Fentanyl), ASTM D5151, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Dental, Healthcare, Laboratory, Veterinarian

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings



\*Typical values have an allowable specification variance  
<sup>†</sup>Contact DASH Medical for a list of tested Chemotherapy drugs