

Nuvo™

White vinyl exam gloves



Super Stretchy, Form Fitting, Opaque White



- Opaque white color
- Super stretch formulation
- Form fitting
- High quality medical grade
- Powder free



"Top Of The Line Product

BEST value glove EVER!! These are very durable but yet flexible for different hand shapes. The size is spot on. The fact the company has fantastic customer service, easy online ordering and the shipments are sent out quickly are just bonuses."

– Starla



The most amazing vinyl exam gloves you'll find. Soft and form fitting, Nuvo™ gloves utilize a unique synthetic vinyl formulation to deliver a comfortable barrier.



PRODUCT DETAILS	
Material	Vinyl
Interior	Powder Free
Grade	Medical Examination
Watertight Test	1.5 AQL
Color	White
Texture	Smooth

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.08 mm 3.1 mil
Finger Thickness	0.06 mm 2.4 mil
Length	230 mm 9.1 in
Tensile Strength (MPA)	11
Elongation %	300

ORDER INFORMATION	
Product Number	NV100
Available Sizes	S, M, L, XL
Glove Count	100/box 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested†	No
Tested Against Fentanyl	No
PFAS Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D5250, STM D5151, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Healthcare, Child Care, Food Service, Veterinarian

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings



800-523-2055 | DASHMEDICAL.COM

DASH_SS_NV100_09.2022

*Typical values have an allowable specification variance
†Contact DASH Medical for a list of tested Chemotherapy drugs