


True Comfort[®] blu

Polychloroprene Exam Gloves



Experience all the feels

Experience the feel of **SILK**



Experience unrivaled **WET GRIP**



Experience a new level of **COMFORT**



true Comfort[®] blu

A new breakthrough in chloroprene polymer technology delivers improved wet grip, increased touch sensitivity and silky smooth comfort.



PRODUCT DETAILS	
Material	Polychloroprene
Interior	Powder Free
Grade	Medical Examination
Watertight Test	1.5 AQL
Color	Ocean Blue
Texture	Fingers

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.06 mm 2.4 mil
Finger Thickness	0.08 mm 3.1 mil
Length	238 mm 9.3 in
Tensile Strength (MPA)	23
Elongation %	500

ORDER INFORMATION	
Product Number	TCB100
Available Sizes	XS, S, M, L, XL
Glove Count	100/box 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested†	Yes
Tested Against Fentanyl‡	No
PFAS Free	Yes
BPA Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D6977, ASTM D6978 (Chemo), ASTM D5151, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Dental, Laboratory, Industrial Safety, Automotive

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings

Experience GIVING BACK.

DASH Medical will make a donation for every case sold of true Comfort blu gloves to organizations that clean and protect the environment. Visit dashtrueblu.com to learn more.



800-523-2055 | DASHMEDICAL.COM

DASH_SS_TCB100_03.2023

*Typical values have an allowable specification variance
†Contact DASH Medical for a list of tested Chemotherapy drugs