

VitalGard®

Latex exam gloves



Top Notch Protection at Rock Bottom Prices



- Aggressively priced
- Excellent range of motion
- Powder free
- Finger textured



"Excellent Buy!"

Excellent product, durable and reliable. And the customer service is the best of any of my suppliers."

- John

VITAL GARD®

LATEX EXAMINATION GLOVES

With so many choices in hand protection these days, maybe you just want to keep it simple: low-cost, high quality. VitalGard® latex gloves are your answer.



PRODUCT DETAILS	
Material	Latex
Interior	Powder Free
Grade	Medical Examination
Watertight Test	2.5 AQL
Color	Natural
Texture	Full

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.11 mm 4.3 mil
Finger Thickness	0.12 mm 4.7 mil
Length	240 mm 9.5 in
Tensile Strength (MPA)	18
Elongation %	650

ORDER INFORMATION	
Product Number	VFP100
Available Sizes	XS, S, M, L, XL
Glove Count	100/box 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested†	No
Tested Against Fentanyl	No
PFAS Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D3578, ASTM D5151, ASTM D5712, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Dental, Healthcare, Laboratory, Veterinarian

PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings



*Typical values have an allowable specification variance
†Contact DASH Medical for a list of tested Chemotherapy drugs