

VitalGard®

Vinyl exam gloves



Top Notch Protection at Rock Bottom Prices



- Clear vinyl exam glove
- Latex free
- Powder free



"We Love Them

We use these gloves at our child care center. We love them! They are reasonably priced, good quality, and they arrive quickly!"
– Savanna

VITAL GARD®

VINYL EXAMINATION GLOVES

Non-latex and powder free solution with a loose comfortable fit for a variety of jobs.



PRODUCT DETAILS

| | |
|-----------------|---------------------|
| Material | Vinyl |
| Interior | Powder Free |
| Grade | Medical Examination |
| Watertight Test | 1.5 AQL |
| Color | Clear |
| Texture | Smooth |

PHYSICAL PROPERTIES (TYPICAL VALUES*)

| | |
|------------------------|-------------------|
| Palm Thickness | 0.08 mm 3.1 mil |
| Finger Thickness | 0.10 mm 3.9 mil |
| Length | 235 mm 9.3 in |
| Tensile Strength (MPA) | 17 |
| Elongation % | 430 |

ORDER INFORMATION

| | |
|-----------------|--------------------------|
| Product Number | VV100 |
| Available Sizes | S, M, L, XL |
| Glove Count | 100/box 10 boxes/case |

CERTIFICATIONS

| | |
|---------------------------|-----|
| Chemotherapy Drug Tested† | No |
| Tested Against Fentanyl | No |
| PFAS Free | Yes |
| FDA Food Contact | Yes |

PRODUCT CONFORMANCE

ASTM D5250, ASTM D5151, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES

Healthcare, Food Service, Janitorial, Veterinarian

| PERFORMANCE LEVEL | LIGHT DUTY | HIGH TACTILITY | MID DUTY | HEAVY DUTY | HIGH RISK |
|-------------------------------|------------------------|-----------------------------|---------------------------------------|------------------------------------|------------------------------|
| CUSTOMER NEED/ APPLICATION | Basic Barrier Coverage | Best Touch Sensitivity | Broad Range of Uses | Increased Durability | Longer Cuff Length |
| | Lowest cost option | Delivers a second skin feel | Perfect for a variety of applications | Great for more rugged applications | Used in higher risk settings |