

ALASTA[®] 100 & 200

Nitrile exam gloves



Extreme Comfort, Precision Fit

Market leader for over 15 years



- Comfortable precision fit, reliable strength
- Available in 100 and 200 glove count boxes
- Chemotherapy drug rated[†]
- Tested against fentanyl



"Love These Gloves!"

I've been using these gloves for 17 years and I love them! They fit my hands perfect."

- Traci

Your Trusted Choice in Gloves - Protecting Hands and People, Made Easy



800-523-2055 | DASHMEDICAL.COM

Alasta[®]

Soft Fit™ Nitrile Gloves

With over 15 years of proven performance, Alasta[®] nitrile exam gloves deliver precision, comfort, and dependable protection powered by Soft Fit™ thin nitrile technology.



PRODUCT DETAILS	
Material	Nitrile
Interior	Powder-Free
Grade	Medical Examination
Watertight Test	1.5 AQL
Color	Violet Blue
Texture	Finger

PHYSICAL PROPERTIES (TYPICAL VALUES*)	
Palm Thickness	0.06 mm 2.4 mil
Finger Thickness	0.08 mm 3.1 mil
Length	245 mm 9.6 in
Tensile Strength (MPA)	35
Elongation %	545

ORDER INFORMATION	
Product Number	ALS100 ALS200
Available Sizes	XS, S, M, L, XL
Glove Count	100/box 10 boxes/case 200/box 10 boxes/case

CERTIFICATIONS	
Chemotherapy Drug Tested†	Yes
Tested Against Fentanyl	Yes
PFAS	Yes
BPA Free	Yes
FDA Food Contact	Yes

PRODUCT CONFORMANCE
ASTM D6319, ASTM D6978 (Chemo and Fentanyl), ASTM D5151, ASTM D6124, FDA 21 CFR 177

RECOMMENDED INDUSTRIES
Dental, Healthcare, Laboratory, Industrial Safety, Automotive, Veterinarian



2 mil

PACKAGED
ALS100: 100 GLOVES/BOX
ALS200: 200 GLOVES/BOX

[Learn More/Request FREE Samples Here](#)



PERFORMANCE LEVEL	LIGHT DUTY	HIGH TACTILITY	MID DUTY	HEAVY DUTY	HIGH RISK
CUSTOMER NEED/ APPLICATION	Basic Barrier Coverage	Best Touch Sensitivity	Broad Range of Uses	Increased Durability	Longer Cuff Length
	Lowest cost option	Delivers a second skin feel	Perfect for a variety of applications	Great for more rugged applications	Used in higher risk settings



800-523-2055 | DASHMEDICAL.COM

DASH_SS_ALS_11.2025

*Typical values have an allowable specification variance
†Contact DASH Medical for a list of tested Chemotherapy drugs